Research Summary – Parkinson's Disease (PD)

Tremor, bradykinesia (i.e., slowness of movement), postural instability and rigidity (i.e., stiffness and resistance to passive movement) are the cardinal motor symptoms of PD.1 In addition, patients also experience non-motor symptoms like dysarthria, hyposmia, sleep disorders, and cognitive impairments. Collectively, these symptoms have a significant impact on the functional ability of patients and their overall quality of life.

Gait analysis and clinical correlations in early Parkinson's disease -

https://pmc.ncbi.nlm.nih.gov/articles/PMC5505527/

- Study set out to identify and quantify spatiotemporal and kinematic gait parameters by 3D gait analysis in a group of PD patients compared to control group.
- Participant size 44 early-stage PD patients (23 male & 21 female) mean age 66.5 ± 9.11 yrs
 & 44 age and sec-matched healthy participants as the control group (22 male & 22 female) mean age 67.00 ± 9.42 yrs.
- Findings
 - \circ Steps per minute 102.46 ± 13.17 steps/min in PD patients vs 113.84 ± 4.30 steps/min in control subjects
 - \circ Stride duration 1.19 \pm 0.18 seconds **right** limb and 1.19 \pm 0.19 seconds **left** limb in PD patients vs 0.426 \pm 0.16 seconds **right** limb and 0.429 \pm 0.23 seconds **left** limb in normal subjects
 - \circ Stance duration 0.74 \pm 0.14 seconds **right** limb and 0.74 \pm 0.16 **left** limb in PD patients vs 13.4 \pm 1.1 seconds **right** limb and 0.83 \pm 0.6 seconds **left** limb in control subjects
 - Swing phase and swing duration differed considerably too (p<0.05), while the stance phase was not statistically significant in patients compared with healthy subjects.
 - \circ Spatial parameters velocity in PD patients (0.082 ± 0.29 m/s) vs healthy subjects (1.33 m/s ± 0.06).
 - Step width, stride length and swing velocity were highly significant parameters, as was average velocity.
 - Reduced step length has been reported as one of the key features of PD gait (Yang et al., 2007)
 - Results demonstrated that early moderate PD patients showed a reduction of spatiotemporal gait parameters (speed and cadence)
 - Loss of balance, short stride and slowness have been found to be the most prevalent clinical features in PD, after exclusively parkinsonian signs (Morris et al., 1994; Morris et al., 1998; Chien et al., 2006; Sofuwa et al., 2005)

Real-World Gait Detection Using a Wrist-Worn Inertial Sensor: Validation Study - https://formative.jmir.org/2024/1/e50035/

- Study set out to validate gait sequence (GS) detection algorithms developed for the wrist position against reference data acquired in a real-world context
- Participant size 83 participants across the health spectrum (PD, MS, heart failure, healthy etc) monitored for 2.5 hours using wrist, lower back and feet sensors. 15 were PD patients.

- In total, 10 algorithms for wrist-based gait detection were validated against a multisensor reference system and compared to gait detection performance using lower back—worn inertial sensors.
- Findings
 - Algorithms applied to the wrist position can detect GSs with high performance in real-world environments.
 - Wrist position is not as accurate as a lower-back position, but significant findings showed a wrist-position still provided accurate information, just at a lower sensitivity.
 - Provides evidence that wrist-worn devices can accurately detect changes in gait.

Development of digital biomarkers for resting tremor and bradykinesia using a wrist-worn wearable device - https://www.nature.com/articles/s41746-019-0217-7

- This study propose a method that sequentially processes epochs of raw sensor data from a single wrist-worn accelerometer by using heuristic and machine learning models in a hierarchical framework to provide continuous monitoring of tremor and bradykinesia
- Findings
 - Results show that sensor derived continuous measures of resting tremor and bradykinesia achieve good to strong agreement with clinical assessment of symptom severity and are able to discriminate between treatment-related changes in motor states.
 - o Approach consists of two steps: Context detection & Symptom severity estimation.
 - The context detection step includes detection of hand movement followed by detection of gait (if hand movement was detected). The symptom severity estimation step includes detection and assessment of resting tremor (if no hand movement was detected) and assessment of bradykinesia (if hand movement was detected and gait was not detected)
 - Hand movement classifier dataset used achieved an accuracy of 92.82% (weighted precision: 93%, recall: 93% and F1 score: 93%) with a false-negative rate of 4.72% and false-positive rate of 10.78%.
 - Gait classifier dataset used achieved an accuracy of 96% (weighted precision: 96%, recall: 96% and F1 score: 96%) with a false-negative rate of 2.45% and false-positive rate of 6.4%.
 - Resting tremor classifier dataset used achieved an accuracy of 83% (weighted precision: 86%, recall: 86% and F1 score: 83%) with a false positive rate of 6.83%

Using wearables to assess bradykinesia and rigidity in patients with Parkinson's disease: a focused, narrative review of the literature -

https://www.researchgate.net/publication/333298943 Using wearables to assess bradykinesia a nd rigidity in patients with Parkinson%27s disease a focused narrative review of the literatur e

- Literature review to consolidate studies that evaluated objective measurement of bradykinesia and rigidity (31 & 8 studies respectively) within PD patients.
- Findings
 - Several studies reported strong associations between wearable-based measures and the gold-standard references for bradykinesia, and, to a lesser extent, rigidity.

TremorSense: Tremor Detection for Parkinson's Disease Using Convolutional Neural Network – https://www.cs.wm.edu/~wsjung/resources/TremorSenseChase.pdf

- Study on TremorSense a PD tremor detection system/application to classify Parkinson's Disease hand tremors.
- Participants 30 PD patients utilizes accelerometers and gyroscopes as wearable sensors on patients' wrists.

TABLE III: Data Segmentation Instances

Label	Rest	Postural	Action	Total
Tremor	15049	1206	4011	20226
No Tremor	17705	504	2017	20226

Findings –

- Self-Evaluation: For self-evaluation, the study trained and tested the model using the data from each patient. Since they trained and tested the model using the same data for each patient, the classification results were supposed to achieve the best score among all evaluation results. As expected, the overall precision, recall, F1 score, and accuracy are 93.31%, 100%, 96.54%, and 96.65%, respectively.
- Cross-evaluation: For cross-evaluation, the study trained and tested the model on using the data from all patient data. They also run 5-fold cross-validation for each round of training and testing. The overall precision, recall, F1 score, and accuracy are 90.64%, 100%, 95.09%, and 95.35%, respectively.
- Leave-one-out-evaluation: For leave-one-out- evaluation, we trained our model using data from 29 patients and tested our model with the data from the single remaining patient. The overall precision, recall, F1 score and accuracy are 88.58%, 100%, 93.95%, and 94.29%, respectively.
- Even though the leave-one-out evaluation accuracy is around 2% worse than selfevaluation, and 1% worse than cross-evaluation, the performance is still strong with an accuracy above 94%.

Parkinson's Disease Action Tremor Detection with Supervised-Leaning Models - https://pmc.ncbi.nlm.nih.gov/articles/PMC10516258/

- The study introduces a PD action tremor detection method to recognize PD tremors from regular activities from TremorSense original dataset.
- Participants 30 PD patients wearing accelerometers and gyroscope sensors on their wrists (18 male & 12 females – mean age of 67.43). The data for the TremorSense dataset was collected using two UG sensor bands containing a three-axis accelerometer and a three-axis gyroscope
- Findings
 - Support Vector Machines (SVMs) perform the best in five-fold cross-validation with over 92% F1 scores. SVMs also show the best performance in the leave-one-out evaluation with over 90% F1 scores.

Classifiers	Accuracy	Precision	Recall	F1-Score
Linear Model	91.37%	91.44%	91.05%	91.24%
KNN	90.21%	90.02%	90.45%	90.51%
CNN	92.24%	92.50%	91.98%	92.21%
SVM	92.28%	92.40%	92.13%	92.26%
TremorSense	83.60%	87.18%	79.27%	83.04%

Datasets – Mobilise-d - https://mobilise-d.eu/data/